

QuickTox™ Kit for DON (Deoxynivalenol) at 0.5 ppm, 1 ppm, or 2 ppm Wheat, Corn & Barley

Highlights:

- Results in 5 minutes
- Simple protocol – no solvents
- Multiple detection levels – one sample prep

Contents of Kit:

- 50 QuickTox Strips packed in a moisture-resistant canister
- 50 reaction vials
- 50 transfer pipettes marked for 1 mL
- Sample Diluent

Items Not Provided:

- Plastic sample cups with lids*
- Graduated cylinder*
- Pipette to dispense 50 μ L*
- Timer

* Available as accessories through EnviroLogix – see list on Page 3



Place 10 grams of sample into cup, add water



Shake vigorously for 30 seconds

Catalog Number AS 104 BG

Intended Use

The QuickTox Kit for DON (deoxynivalenol, vomitoxin) is designed to quickly extract and screen wheat, corn and barley for the presence of deoxynivalenol residues. This Kit will provide a qualitative screen for deoxynivalenol residues at any of three cutoff levels in bulk grain samples: 2 ppm 1 ppm and 0.5 ppm.

How the Test Works

A composite grain sample is first collected, then extracted to solubilize any deoxynivalenol present. Each sample should be ground to a fineness of at least 20 mesh and extracted with tap water. This extract is further diluted in the Sample Diluent provided for testing with the QuickTox Kit.

Each QuickTox Strip has an absorbent pad at each end. The protective tape with the arrow indicates which end of the strip to insert into the reaction vial. The sample extract travels up the membrane strip and is absorbed into the larger pad at the top of the strip. The portion of the strip between the protective tape and the absorbent pad at the top of the strip is used to view the reactions as described under “Interpreting the Results.”

Preparation of the Sample

Step 1: Determine Number and Size of Sub-samples

1. Collect a composite sample according to your own sampling plan or USDA/GIPSA guidelines. Consult USDA/GIPSA reference documents such as <http://archive.gipsa.usda.gov/reference-library/handbooks/grain-insp/grbook1/bk1.pdf> to help design a plan that fits your needs.
2. Grind samples using a mill which provides a sample that passes through at least a 20 mesh sieve. Mix ground material thoroughly before sub-sampling.

Step 2: Extract sample

1. Weigh 10 grams of milled sample into a disposable sample cup with lid and add five volumes of room temperature tap water (e.g. for 10 grams, add 50 mL).
2. Cap sample cup tightly and shake vigorously by hand for 30 seconds. Ensure that entire sample appears completely wet and thoroughly mixed.
3. Allow sample to settle until 2 distinct layers are visible and fine particles are mostly settled (1-3 minutes, depending on grind). The top layer containing the DON residues will be used in testing.




Allow extract to settle into 2 layers, use top layer and avoid particles



Add Sample Diluent to vial, then add extract and mix



Place strip in vial



Step 3: Dilute extract with Sample Diluent (select testing level)

	0.5 ppm	1.0 ppm	2.0 ppm
1.	Using a transfer pipette (provided), place 1 mL Sample Diluent into reaction vial.	Using a transfer pipette (provided), place 1 mL Sample Diluent into reaction vial.	Using a transfer pipette (provided), place 1 mL Sample Diluent into reaction vial.
See "Notes" section (p. 4) for directions on using the transfer pipette to measure correctly.			
2.	Using a pipette and clean tip, transfer 200 µL from the top layer of sample into reaction vial containing Sample Diluent. Avoid particulates.	Using a pipette and clean tip, transfer 100 µL from the top layer of sample into reaction vial containing Sample Diluent. Avoid particulates.	Using a pipette and clean tip, transfer 50 µL from the top layer of sample into reaction vial containing Sample Diluent. Avoid particulates.
3.	Mix diluent and sample extract thoroughly by stirring with the tip of the pipette	Mix diluent and sample extract thoroughly by stirring with the tip of the pipette	Mix diluent and sample extract thoroughly by stirring with the tip of the pipette

NOTE: Pipette with care to avoid contamination and ensure correct volumes are used to prepare the test samples. Do not reuse diluted samples. Use a new disposable transfer pipette, pipette tip, and reaction vial for each sample.


How to Run the QuickTox Strip Test

Note: Allow QuickTox Strips and Sample Diluent to warm up to room temperature before preparing samples.

1. Allow refrigerated canisters to come to room temperature before opening. Remove the QuickTox Strips to be used. Avoid bending the strips. Reseal the canister immediately.
2. Place the strip into the reaction vial containing the diluted sample extract. The arrow tape on the end of the strip should point into the reaction vial.
3. The sample extract will travel up the strip. Reaction vials will stand on their own or may be inserted into the cardboard rack provided.
4. Allow the strip to develop for 5 minutes before making final assay interpretations. Strips should be read **promptly at 5 minutes** while wet.
5. To retain the strip, cut off the strip immediately below the bottom line (as indicated below). Discard the bottom section of the strip covered by the arrow tape.

Interpreting the Results

The QuickTox Strip for DON has three lines. The top line is a Negative Reference (NR) that develops a signal similar to a Negative DON sample. The bottom line is a Test Reference Line (REF) that develops a signal similar to how the Test Line should appear at the desired testing level. The middle line is the Test Line and is compared to the bottom REF Line. Results are invalid if either Reference Line fails to develop.




If the middle Test Line color intensity is:	Screening at 0.5 ppm	Screening at 1.0 ppm	Screening at 2.0 ppm
	The amount of DON present in the sample is:		
Darker than the (Bottom) REF Line	less than 0.5 ppm	less than 1 ppm	less than 2 ppm
Equal to or lighter than the (Bottom) REF Line	0.5 ppm or greater	1 ppm or greater	2 ppm or greater

Kit Storage

This QuickTox Kit should be stored refrigerated. Note the shelf life on the kit box. The kit may be used in field applications; however, prolonged exposure to high temperatures may adversely affect the test results. Do not open the desiccated canister until ready to use the strips.

Precautions and Notes

- This kit is designed to screen for presence or absence only, and is not designed to be quantitative.
- This product is currently not applicable for use in testing crops other than those listed in the Intended Use section.
- As with all screening tests, it is recommended that results be confirmed by an alternate method when necessary.
- The assay has been optimized for use with the protocol provided in the kit. Deviation from this protocol may invalidate the results of the test.
- The results generated through the proper use of this diagnostic tool reflect the condition of the working sample directly tested. Extrapolation as to the condition of the originating lot, from which the working sample was derived, should be based on sound sampling procedures and statistical calculations which address random sampling effects, non-random lot sampling effects and assay system uncertainty. A negative result obtained when properly testing the working sample does not necessarily mean the originating lot is entirely negative for the analyte in question.
- Results should be read promptly at 5 minutes. Results read beyond this time may be less reliable.
- Protect all components from hot or cold extremes of temperature when not in use. Do not leave in direct sunlight or in vehicle.



- To use disposable transfer pipettes:
 - Squeeze bulb tightly and insert tip in the Sample Diluent
 - Release pressure to draw liquid up past the 1 mL mark
 - Squeeze carefully to expel excess Diluent back into the bottle so that the liquid left in the pipette is at the 1 mL mark
 - Move the pipette over to the reaction vial and expel the 1 mL Diluent.

- For convenience, some accessories can be ordered through EnviroLogix (see list below).

Optional Items Available:

- Graduated cylinder ACC 023
- Set of 50 sample cups with caps ACC 012
- 50 μ L fixed volume Mini-Pet ACC 051
- 50 μ L pipette tips ACC 053





**For Technical Support
Contact Us At:**

EnviroLogix
500 Riverside Industrial
Parkway
Portland, ME 04103-1486
USA

Tel: (207) 797-0300
Toll Free: 866-408-4597
Fax: (207) 797-7533

e-mail:
info@envirologix.com

website:
www.envirologix.com



LIMITED WARRANTY

EnviroLogix Inc. ("EnviroLogix") warrants the products sold hereunder ("the Products") against defects in materials and workmanship when used in accordance with the applicable instructions for a period not to extend beyond a product's printed expiration date. If the Products do not conform to this Limited Warranty and the customer notifies EnviroLogix in writing of such defects during the warranty period, including an offer by the customer to return the Products to EnviroLogix for evaluation, EnviroLogix will repair or replace, at its option, any product or part thereof that proves defective in materials or workmanship within the warranty period.

ENVIROLOGIX MAKES NO OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. The warranty provided herein and the data, specifications and descriptions of EnviroLogix products appearing in EnviroLogix published catalogues and product literature are EnviroLogix' sole representations concerning the Products and warranty. No other statements or representations, written or oral, by EnviroLogix' employees, agents or representatives, except written statements signed by a duly authorized officer of EnviroLogix Inc., are authorized; they should not be relied upon by the customer and are not a part of the contract of sale or of this warranty.

EnviroLogix does not warrant against damages or defects arising in shipping or handling, or out of accident or improper or abnormal use of the Products; against defects in products or components not manufactured by EnviroLogix, or against damages resulting from such non-EnviroLogix made products or components. EnviroLogix passes on to customer the warranty it received (if any) from the maker thereof of such non-EnviroLogix made products or components. This warranty also does not apply to Products to which changes or modifications have been made or attempted by persons other than pursuant to written authorization by EnviroLogix.

THIS WARRANTY IS EXCLUSIVE. The sole and exclusive obligation of EnviroLogix shall be to repair or replace the defective Products in the manner and for the period provided above. EnviroLogix shall not have any other obligation with respect to the Products or any part thereof, whether based on contract, tort, strict liability or otherwise. Under no circumstances, whether based on this Limited Warranty or otherwise, shall EnviroLogix be liable for incidental, special, or consequential damages.

This Limited Warranty states the entire obligation of EnviroLogix with respect to the Products. If any part of this Limited Warranty is determined to be void or illegal, the remainder shall remain in full force and effect.

License

EnviroLogix has developed this kit using proprietary reagents.

EnviroLogix, the EnviroLogix logo, and QuickTox are trademarks of EnviroLogix Inc.

Patent Pending

© EnviroLogix 2008