

Highlights:

- Results in 3 minutes
- Quantitative and traceable results in QuickScan
- Quantitate up to 12 ppm with no extra steps
- Range extension in corn and wheat up to 30 ppm with additional dilution step.
- Read strips wet – no drying necessary
- Simple protocol
- No incubation equipment needed

Contents of Kit:

- 50 QuickTox Strips packed in a moisture-resistant canister
- 100 reaction vials
- 100 pipette tips
- DB6 Buffer, kit lot specific

Items Not Provided:

- Orbital/rotary shaker
- Plastic sample cups with lids* or other extraction vessels
- Bottled, distilled or deionized water
- 20 mesh screen
- Graduated cylinder*
- Pipettes and tips to deliver 200 µL and 800 µL*
- Approved Coffee Filters (EnviroLogix validated)*
- Timer
- Scissors
- Osterizer Blender Model 4094 or equivalent (wheat gluten only)
- Oster Mini Blend Jar Part no. 4888-3 (for samples <35g, wheat gluten only)
- Microcentrifuge (wheat gluten only)*
- Tubes and pipettes for centrifugation (wheat gluten only)*
- QuickScan System*

*Available as accessories – see list on Page 7

Catalog Number AQ 254 BG

Part #11176

Intended Use

The QuickTox Kit for QuickScan DON3 is designed to quickly extract and screen milled samples for the presence of Deoxynivalenol (DON) residues. The QuickTox Kit will then provide quantitative results when used in conjunction with the QuickScan System. The QuickTox Kit is designed to provide quantitative results up to 12 ppm in the base range and up to 30 ppm in corn and wheat with additional sample dilution. Commodities include corn, corn flour, corn germ, corn gluten meal, corn gluten feed, DDGS; wheat, wheat bran, wheat midds, wheat flour, white wheat flour, wheat gluten, wheat red dog; barley, malted barley; oats; sorghum; soybean meal; milled rice and rough rice; and whole rye.

How the Test Works

A composite sample is first collected, then extracted to solubilize any DON present. Each sample should be ground to a fineness of 20 mesh and extracted with room temperature distilled, deionized, flat bottled water, or qualified drinkable (potable) tap water. This extract is further diluted for testing with the QuickTox Kit.

Each QuickTox Strip has an absorbent pad at each end. The protective tape with the arrow indicates which end of the strip to insert into the reaction vial. The sample extract travels up the membrane strip and is absorbed into the larger pad at the top of the strip. At three minutes, the strip is cut off at the top of the arrow tape, the bottom pads are discarded, and the strip is inserted into the QuickScan reader to obtain quantitative results.



Preparation of the Sample

Please note: sample extract should be tested immediately after dilution with DB6 Buffer (Step 7). Make sure strips and DB6 Buffer are at room temperature and ready for use before the dilution step.

Determine number and size of sub-samples

1. Collect a composite sample according to your own sampling plan or USDA/GIPSA guidelines. Consult USDA/GIPSA reference documents such as www.gipsa.usda.gov/fgis/handbook/gihbk1_inspec.aspx to help design a plan that fits your needs.
2. Grind samples using a mill which provides a sample such that ≥95% passes through a 20 mesh sieve. Mix ground material thoroughly before sub-sampling.



Correct 20 mesh grind




Measure water, add to ground sample (Step 3)



Shake mechanically or by hand (Step 4)



Filter extract (Step 5)

Alternate Step 5 (wheat+corn only)

Allow to settle into layers; angling the cup to draw off extract will help avoid the top foamy layer and the lower particulate layer

Extract Sample

Different commodities require different extraction ratios. Choose your commodity and follow the instructions for extraction listed in the table below:

Note: Potable tap water may be utilized instead of distilled, deionized or bottled water, following customer validation of water supply. Please contact EnviroLogix to purchase a control set and protocol that can be used to verify your water supply.

Corn gluten meal, corn gluten feed, wheat flour, malted barley, DDGS, whole rye	Corn, corn flour, corn germ, wheat, wheat midds, wheat bran, wheat red dog, barley, rough rice, soybean meal, wheat gluten†	Oats*, white wheat flour*, milled rice, sorghum
4X	5X	6X
3. Weigh 20 to 50 grams of milled sample into a disposable sample cup with lid or other suitable container and add four volumes of room temperature bottled (flat), deionized, or distilled water (4 mL per gram of sample, i.e. 20 grams, add 80 mL).	3. Weigh 20 to 50 grams of milled sample into a disposable sample cup with lid or other suitable container and add five volumes of room temperature bottled (flat), deionized, or distilled water (5 mL per gram of sample, i.e. 20 grams, add 100 mL).	3. Weigh 20 to 50 grams of milled sample into a disposable sample cup with lid or other suitable container and add six volumes of room temperature bottled (flat), deionized, or distilled water (6 mL per gram of sample, i.e. 20 grams, add 120 mL).

* Shake time one minute; see Step 4

† Blender extraction + centrifugation required – see Blender Extraction instructions below

- Cap sample cup tightly, disperse water making sure the entire sample is wetted and place on shaker (alternatively, shake vigorously by hand) for 30 seconds (*oats and white wheat flour requires a 1 minute shake). Shaker should be operated at the highest speed. Samples that are not thoroughly mixed may adversely affect test results due to incomplete extraction.
- Add an approved coffee filter (example: BUNN part #BUNBCF100B) to a clean vessel and pour extract into the filter. Pull back an edge of the filter to gain access to the filtered extract.

Alternatively, for corn and wheat only, allow extract to separate into 2 layers. The top layer containing the DON residues will be used in testing. In some instances a foamy layer will float above the desired top layer. The best technique to retrieve this extract for Step 7 below is to tip the extraction cup at a 45 degree angle, exposing the supernatant beneath the foamy layer, avoiding particulates.

Blender Extraction + centrifugation (wheat gluten):

- Determine weight of sample and calculate appropriate volume of water (5 mL per gram of sample, i.e. 50 grams, add 250 mL). Choose blender vessel based upon sample weight:
 - For samples between 20 and 35g, use the Oster mini-blend jar part number 4888-3.
 - For samples between 35 and 50g, use the standard glass extraction vessel that arrived with the blender.
- Add the appropriate volume of water to the blender vessel first, then add the sample.

Blender extraction (wheat gluten only)

Add water and wheat gluten to blender and blend for 1 minute



Immediately centrifuge for 1 minute and use clear liquid top layer for testing

Dispense liquids:

Add Buffer to Dilution Vial first, then add extract; mix well with pipette tip (Step 6,7)



Transfer mixture to Testing Vial (Step 8)



Insert strip
Wait 3 minutes for results

Cut strip and place in QuickScan reader immediately—no drying step!

- Blend on highest speed for 1 minute.
- Immediately fill a microcentrifuge tube with extract. Centrifuge for 1 minute at 2000 x g (not RPM). The clear liquid top layer containing the DON residues will be used in testing.

Dispense liquids (use 2 vials and 2 pipette tips for each sample being tested: the first set for dilution, the second set for testing):

Dilution (vial #1):

- Using a calibrated pipette with a **new tip**, place 800 microliters (800 μ L) DB6 Buffer into the (first) dilution vial. Take care not to contaminate DB6 Buffer—use a new tip for each test and keep covered when not in use. **Please note:** DB6 Buffer is matched with specific DON3 kit lot numbers; be sure to use the DB6 that is provided with the kit (do not mix and match buffers with different kit lots). There is a "use with" label on the DB6 that will indicate the matching DON3 lot number.
- With **another new** pipette tip, remove 200 μ L of the extract. Add that extract to the dilution vial containing the DB6 Buffer and mix well with pipette by stirring or drawing liquids up and down in the pipette tip until the mixture is uniform.

Testing (vial #2):

- With the same pipette tip, add 200 μ L of the mixed sample to the (second) reaction vial, used for testing.

NOTE: Samples that are not thoroughly mixed and/or accurately pipetted will adversely affect test results. When following the alternate settling protocol for corn and wheat, **avoid foam and particulates during pipetting, and ensure that the pipette tip does not become clogged with particulate.** The run volume in the reaction vial should be 200 μ L. Do not reuse diluted samples. Use new reaction vials for each sample. Use two pipette tips (one for DB6 Buffer, one for extract) for each sample.


For testing corn and wheat samples at levels greater than 12 ppm (up to 30 ppm):

If after running and reading the test, the initial result is greater than 12 ppm ("> 12 ppm" on QuickScan), and further knowledge about the level of contamination is desired, samples can be retested by further dilution of the sample extract.

- In a separate dilution vial (not provided), combine seven parts of water with one part filtered extract (for example, 700 μ L water and 100 μ L filtered DON extract). **Mix well.**
- Using a calibrated pipette with a **new tip**, place 800 μ L DB6 Buffer into a dilution vial.
- With **another new** pipette tip, add 200 μ L of the newly diluted extract to the dilution vial containing Buffer and mix thoroughly.
- With the same pipette tip, add 200 μ L of the mixed sample to a new reaction vial, used for testing. Follow the instructions under How to Run the QuickTox Strip Test. Choose 1:7.5 under the dilution tab on the QuickScan Results Screen—the System will calculate and record the DON level in the diluted samples.

How to Run the QuickTox Strip Test

- Allow refrigerated canisters to come to room temperature before opening. Remove the QuickTox Strips to be used. Avoid bending the strips. Reseal the canister immediately.



- Place the strip into the reaction vial containing the 200 μ L of DB6 buffer/sample mixture (Vial #2). The arrow tape on the end of the strip should point into the reaction vial.
- The sample extract will travel up the strip (flow may not be visible immediately, this is normal). Reaction vials will stand on their own.
- Allow the strip to develop for **three (3)** minutes. Immediately cut off and discard the bottom section of the strip covered by the arrow tape and insert strip into the QuickScan reader for quantitation.

Use of the QuickScan System

Detailed instructions for use of the QuickScan System are supplied with each unit, and can also be found at envirologix.com/quickscan. In summary, a strip is inserted into the reader and the strips are read by touching or clicking on the “Read Test” area of the screen. Results are then recorded in an electronic worksheet, allowing each user to report and track data easily.

Results are reported in the range of 0 to 12 ppm. The limit of detection (LOD) is 0.29 ppm; accuracy of results less than 0.29 ppm should not be assumed. Results greater than 12 ppm are reported as "> 12 ppm."

Kit Storage

This QuickTox Kit should be stored refrigerated. Note the shelf life on the kit box. Prolonged exposure to high temperatures may adversely affect the test results. Do not open the desiccated canister until ready to use the strips.

Cross-reactivity

The following mycotoxins have been tested with this kit and no false positive results occurred at the 100 ppm level: Aflatoxin B₁, Fumonisin B₁, Ochratoxin A, and Zearalenone.


Precautions and Notes

- Running for 3 minutes and reading the strip(s) promptly is required for accurate results.
- Accuracy of results less than 0.29 ppm should not be assumed.
- This product is currently not applicable for use in testing any other grains.
- Pipettes lose calibration accuracy over time. Calibrate or replace pipettes at least annually.
- Observe any applicable regulations when disposing of samples and extracts.
- This assay is calibrated against wheat and corn reference samples supplied by Trilogy Analytical Laboratory, Washington, MO, and other vendors and associated HPLC data. Where possible, performance in other sample matrices has been validated using naturally contaminated samples. Where naturally contaminated samples are not available, performance has been validated using fortified samples.
- As with all screening tests, it is recommended that results be confirmed by an alternative method when necessary.
- The assay has been optimized for use with the protocols provided in the kit. Deviation from these protocols may invalidate the results of the test. Room-temperature components, proper and thorough mixing, accurate pipetting, and using the correct corresponding DB6 Buffer provided in the kit are essential to



*Place strip in
QuickScan carrier*





accurate results.

- The results generated through the proper use of this diagnostic tool reflect the condition of the working sample directly tested. Extrapolation as to the condition of the originating lot, from which the working sample was derived, should be based on sound sampling procedures and statistical calculations which address random sampling effects, non-random seed lot sampling effects and assay system uncertainty. A negative result obtained when properly testing the working sample does not necessarily mean the originating lot is entirely negative for the analyte in question.
- Protect all components from hot or cold extremes of temperature when not in use. Do not leave in direct sunlight or in vehicle.
- For convenience, accessories can be ordered from EnviroLogix (see list, next page).

Accessories:

Available through EnviroLogix:	Catalog No.	Part #
▪ QuickScan™ System	ACC 331	12721
▪ Sample cups with lids (500/case) <i>for samples up to 30 g; larger samples require different mixing vessels</i>	ACC 012-CS	10167
▪ Graduated cylinder (100 mL)	ACC 068	11207
▪ 200 µL MiniPet pipette (one/location free)	ACC 067	11206
▪ Coffee filters (100)	ACC 083	11434
▪ Microcentrifuge	ACC 064E	11204
▪ Centrifugation Set: Disposables for 50 tests	ACC 010	11214
▪ QuickTox Dilution Set <i>(enough for 100 dilutions)</i>	ACC 080	11219
▪ 1 mL adjustable pipette <i>(helpful for extended dilution)</i>	ACC 1303	11964
▪ Pipette tips for 1 mL pipette (50)	20-0107	12243





**For Technical Support
Contact Us At:**

EnviroLogix
500 Riverside Industrial
Parkway
Portland, ME 04103-1486
USA
Tel: (207) 797-0300
Toll Free: 866-408-4597
Fax: (207) 797-7533

e-mail:
info@envirologix.com

website:
www.envirologix.com



LIMITED WARRANTY

EnviroLogix Inc. ("EnviroLogix") warrants the products sold hereunder ("the Products") against defects in materials and workmanship when used in accordance with the applicable instructions for a period not to extend beyond a product's printed expiration date. If the Products do not conform to this Limited Warranty and the customer notifies EnviroLogix in writing of such defects during the warranty period, including an offer by the customer to return the Products to EnviroLogix for evaluation, EnviroLogix will repair or replace, at its option, any product or part thereof that proves defective in materials or workmanship within the warranty period.

ENVIROLOGIX MAKES NO OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING BUT NOT LIMITED TO ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. The warranty provided herein and the data, specifications and descriptions of EnviroLogix products appearing in EnviroLogix published catalogues and product literature are EnviroLogix' sole representations concerning the Products and warranty. No other statements or representations, written or oral, by EnviroLogix' employees, agents or representatives, except written statements signed by a duly authorized officer of EnviroLogix Inc., are authorized; they should not be relied upon by the customer and are not a part of the contract of sale or of this warranty.

EnviroLogix does not warrant against damages or defects arising in shipping or handling, or out of accident or improper or abnormal use of the Products; against defects in products or components not manufactured by EnviroLogix, or against damages resulting from such non-EnviroLogix made products or components. EnviroLogix passes on to customer the warranty it received (if any) from the maker thereof of such non-EnviroLogix made products or components. This warranty also does not apply to Products to which changes or modifications have been made or attempted by persons other than pursuant to written authorization by EnviroLogix.

THIS WARRANTY IS EXCLUSIVE. The sole and exclusive obligation of EnviroLogix shall be to repair or replace the defective Products in the manner and for the period provided above. EnviroLogix shall not have any other obligation with respect to the Products or any part thereof, whether based on contract, tort, strict liability or otherwise. Under no circumstances, whether based on this Limited Warranty or otherwise, shall EnviroLogix be liable for incidental, special, or consequential damages.

This Limited Warranty states the entire obligation of EnviroLogix with respect to the Products. If any part of this Limited Warranty is determined to be void or illegal, the remainder shall remain in full force and effect.

License

EnviroLogix has developed this kit using proprietary reagents.

EnviroLogix, the EnviroLogix logo, QuickTox, and QuickScan are trademarks of EnviroLogix Inc.

© EnviroLogix 2019



Revision nr.2
Dated 01/01/2019
Page n. 1 / 5

Safety data sheet

SECTION 1. Identification of the substance/mixture and of the company/undertaking

1.1 Product identifier	DB 6 Dilution Buffer
Trade name:	11151 (KR-268)
Part number:	
1.2 Relevant identified uses of the substance or mixture and uses advised against application of the substance / the preparation :	Laboratory chemicals; kit component. Not to be used for purposes other than those specified in product literature.
1.3 Details of the supplier of the safety data sheet	
Manufacturer/Supplier:	EnviroLogix Inc., 500 Riverside Industrial Pkwy, Portland ME 04103, USA Phone: (207) 797-0300
1.4 Emergency telephone number:	(207) 797-0300 Technical Service

SECTION 2. Hazards identification.

2.1 Classification of the substance or mixture	
Classification according to 29CFR 1910.1200:	Not Classified
2.2 Label elements	
Labeling according to 29CFR 1910.1200	
Pictogram:	None
Signal word:	None
Hazard Statements:	None
2.3 Other Statements:	None

SECTION 3. Composition/information on ingredients.

Chemical name	CAS No	EC No	Classification According to 29CFR 1910.1200	Amount (%)
Sodium Tetraborate Decahydrate	1303-96-4	215-540-4	H360 Rep 1B	< 3 %

SECTION 4. First aid measures.

4.1 Description of first aid measures	
After inhalation :	<u>In case of inhalation:</u> Remove to fresh air. If not breathing give artificial respiration. Get medical attention immediately.
After skin contact :	<u>In case of skin contact:</u> Remove contaminated clothing and shoes immediately. Wash affected area with mild soap or detergent for at least 10 minutes or until no evidence of chemical remains.
After eye contact :	<u>In case of eye contact:</u> immediately flush eyes with plenty of water for at least 15 minutes. Lifting eyelids occasionally, until no evidence of chemical remains. Get medical attention immediately.
After swallowing :	<u>In case of ingestion:</u> DO NOT induce vomiting unless directed to do so by medical personnel. Never give anything by mouth to an unconscious person. Call a physician immediately.
4.2 Most important symptoms and effects, both acute and delayed:	None
4.3 Indication of any immediate medical attention and special treatment needed:	None

SDS DB6 Dilution Buffer



Revision nr.2
Dated 01/01/2019
Page n. 2 / 5

SECTION 5. Firefighting measures.

5.1 Extinguishing media:	CO2, extinguishing powder or water spray. Fight larger fires with water spray or alcohol resistant foam.
5.2 Special hazards arising from the substance or mixture:	None
5.3 Advice for firefighters:	Wear protective gear appropriate for fire conditions including respiratory protective gear.

SECTION 6. Accidental release measures.

6.1 Personal precautions, protective equipment and emergency procedures:	In the case of spilled mixture wear gloves to prevent skin contact. In the case of a large spill, additional protection is recommended.
6.2 Environmental precautions:	Do not discharge mixture to sewer system or waterways.
6.3 Methods and material for containment and cleanup:	Absorb in paper towel or suitable absorbent for larger spills and discard in appropriate waste. Clean with water afterwards.
6.4 References to other sections:	For safe handling refer to Section 7. For information on PPE refer to Section 8. For disposal refer to Section 13

SECTION 7. Handling and storage.

7.1 Precautions for safe handling:	Practice good chemical hygiene when handling. Avoid contact with eyes, skin, and clothing.
7.2 Conditions for safe storage, including any incompatibilities:	Store in tightly closed, non-metal container, in a corrosive compatible area. Prevent direct sunlight and heat. Store in well aired storage rooms.
7.3 Specific end use(s):	Apart from the uses mentioned in section 1.2, no other specific uses are stipulated

SDS DB6 Dilution Buffer



Revision nr.2
Dated 01/01/2019
Page n. 3 / 5

SECTION 8. Exposure controls/personal protection.

8.1 Exposure limits:	Components with limit values that require monitoring at the workplace:	EH40/2005 8 Hr TWA = 5mg/m ³	OSHA 8 Hr TWA = 10 mg/m ³
8.2 Exposure Controls:			
8.2.1 Engineering controls	Facilities using this mixture should be equipped with an eyewash and safety shower. Use general or local exhaust ventilation to keep airborne concentrations below permissible exposure limits.		
8.2.2 General protective and hygienic measures:	The usual precautionary measures should be adhered to when handling chemicals.		
Eye Protection:	Safety glasses with side shields, goggles. Use equipment for eye protection tested and approved under appropriate government standards such as NIOSH (US) or EN 166 (EU). Eye and face protection regulations are described by OSHA (US) in 29CFR1910.133. Do not wear contact lenses when working with chemicals		
Hand Protection:	Handle with gloves. Gloves must be inspected prior to use. Use proper glove removal technique (without touching glove's outer surface) to avoid skin contact with this product. Dispose of contaminated gloves after use in accordance with applicable laws and good laboratory practices. Wash and dry hands. The selected protective gloves have to satisfy the specifications of EU Directive 89/686/EEC and the standard EN 374 derived from it.		
Breathing Equipment:	Appropriate respiratory protection should be determined according to local conditions using risk analysis protocols. An approved disposable air purifying particulate respirator may be used as a backup to engineering controls. Always use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU).		
8.2.3 Environmental exposure controls:	Contain spills, do not allow into environment		

SECTION 9. Physical and chemical properties.

9.1 Information on basic physical and chemical properties:	
a) Appearance:	Clear liquid, colorless to slight yellow.
b) Odor:	None
c) Odor Threshold:	No Data Available
d) pH:	8.6
e) Melting point/freezing point:	No Data Available
f) Boiling point/Boiling range:	No Data Available.
g) Flash point:	Not applicable.
h) Evaporation rate:	No Data Available
i) Flammability (solid, gaseous):	No Data Available
j) Upper/lower flammability or explosive limits:	No Data Available
k) Vapor pressure:	No Data Available
l) Vapor density:	No Data Available
m) Relative density:	No Data Available
n) Solubility(ies):	Fully miscible, water.
o) Partition Coefficient: n-Octanol/water:	No Data Available
p) Auto-ignition temperature:	No Data Available
q) Decomposition temperature:	No Data Available
r) Viscosity:	No Data Available
a) Explosive properties:	No Data Available.
t) Oxidizing properties:	No Data Available
9.2 Other information:	No further relevant information available.

SDS DB6 Dilution Buffer



Revision nr.2
Dated 01/01/2019
Page n. 4 / 5

SECTION 10. Stability and reactivity.

10.1 Reactivity:	No data available
10.2 Chemical Stability:	Stable under normal temperatures and pressures.
10.3 Possibility of hazardous reactions:	Under normal conditions of storage and use, hazardous reactions will not occur.
10.4 Conditions to avoid:	No specific data
10.5 Incompatible materials:	No Data Available.
10.6 Hazardous decomposition products:	Under normal conditions of storage and use, hazardous decompositions products should not be produced.

SECTION 11. Toxicological information.

Information on Toxicological Effects	
Acute effects (toxicity tests):	No Data Available
Sensitization:	No sensitizing effects known
CMR (carcinogenicity, mutagenicity and toxicity for reproduction) effects:	No CMR effects.
Additional toxicological information:	No Additional Information
SECTION 12. Ecological information.	
12.1 Toxicity:	No Data Available
12.2 Persistence and degradability :	No Data Available
12.3 Bio accumulative potential:	No Data Available
12.4 Mobility in soil :	No Data Available
12.5 Results of PBT and vPvB assessment:	Not available as a chemical safety assessment, not required/not conducted.
12.6 Other adverse effects:	No Data Available

SECTION 13. Disposal considerations.

Waste treatment methods:	Contact a licensed professional waste disposal service to dispose of this material. Disposal of surplus or waste solutions must be in accordance with applicable local, state, and national laws and regulations.
--------------------------	---

SECTION 14. Transport information.

SDS DB6 Dilution Buffer



Revision nr.2
Dated 01/01/2019
Page n. 5 / 5

14.1 UN-Number DOT, ADR, ADN, IMDG, IATA : Not Hazardous for Transport
 14.2 UN proper shipping name DOT, ADR, ADN, IMDG, IATA : Not Hazardous for Transport
 14.3 Transport hazard classes DOT, ADR, ADN, IMDG, IATA: Not Hazardous for Transport
 14.4 Packing group (DOT, ADR, IMDG, IATA): Not Hazardous for Transport
 14.5 Environmental hazards No environmental hazard.
 14.6 Special precautions for user : None
 14.7 Transport in bulk according to Annex II of MARPOL73/78 and the IBC code: No information available.

SECTION 15. Regulatory information.

15.1 Safety, health, and environmental Regulations Not a hazardous material
 US Federal Regulations
 OSHA Not listed
 SARA 313
 US State Regulations Not Listed
 European/International Regulations
 European labeling in accordance with EC Directives Not hazardous according to European directives
 15.2 Chemical Safety Assessment Not carried out

SECTION 16. Other information.

This information is true based on our present knowledge. However, EnviroLogix makes no representation of its accuracy or completeness. Persons receiving this information must exercise their independent judgment in determining the product's safety and suitability for its intended use. This document shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship

EHS Department
EnviroLogix Inc.

Codes:

H360 May damage fertility or the unborn child

SDS DR6 Dilution Buffer

