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**ENVIROLOGIX**

**Safety data sheet**

**SECTION 1. Identification of the substance/mixture and of the company/undertaking**

**1.1 Product identifier**  
Trade name: Extraction Buffer  
Part number: EB17(1198, 12382)

**1.2 Relevant identified uses of the substance or mixture and uses advised against:**  
Laboratory chemicals; kit component. Not to be used for purposes other than those specified in product literature.

**1.3 Details of the supplier of the safety data sheet**  
Manufacturer/Supplier: EnviroLogix Inc., 500 Riverside Industrial Pkwy, Portland ME 04103, USA  
Information department: Technical Service  
1.4 Emergency telephone number: (207) 797-0300

**SECTION 2. Hazards identification.**

**2.1 Classification of the substance or mixture**  
Classification according to OSHA 29CFR 1910.1200 and Regulation EC 1272/2008 (CLP):

Flammable Solid category 2	H228	Flammable solid
Acute Toxicity Oral 4	H302	Harmful if swallowed
Acute Toxicity Inhalation 4	H322	Harmful if inhaled
Skin Irritation category 2	H315	Causes skin irritation
Serious eye damage category 1	H318	Causes serious eye damage
Specific Target Organ Toxicity Single Exposure category 3	H335	May cause respiratory irritation
Aquatic Toxicity-Chronic category 3	H412	Harmful to the environment with long lasting effects

**2.2 Label elements**  
Labeling according to OSHA 29CFR 1910.1200 and Regulation (EC) 1272/2008

Hazard pictograms:

Signal word: Danger

Hazard statements:  
H228 Flammable solid.  
H302 + H322 Harmful if swallowed or inhaled  
H315 Causes skin irritation.  
H318 Causes serious eye damage.  
H335 May cause respiratory irritation.  
H412 Harmful to aquatic life with long lasting effects.

Precautionary statements:  
P264 Wash hands thoroughly after handling.  
P273 Avoid release to the environment.  
P280 Wear protective gloves/ eye protection.  
P301 + P312 IF SWALLOWED: Call a POISON CENTER/doctor/physician if you feel unwell.  
P304 + P340 IF INHALED: Remove to fresh air and keep comfortable for breathing.  
P305 + P351 + P338 If in Eyes: Rinse cautiously with water for several minutes; remove contact lenses if present and easy to do. Continue rinsing.  
P403 + P233 Store in a well ventilated place. Keep container tightly closed.

**2.3 Other hazards:** No additional hazards listed

SDS: EB17

Revision nr.2  
Dated 04/21/2019  
Page n. 2 / 6

**ENVIROLOGIX**

**SECTION 3. Composition/information on ingredients.**

**3.1 Substances:** Information not relevant

**3.2 Mixtures:** Extraction Reagent Powder (EB17)

Chemical name	CAS No	EC No	Amount (%)	Classification OSHA 29CFR1910.1200
Sodium Lauryl Sulfate	151-21-3	205-788-1	60 to 85	Flam. Sol. 2 H228; Acute Tox. Oral 4 H302; Acute Tox. Inhal. 4 H322; Skin Irrit. 2 H315; Eye Dam. 1 H318; STOT SE 3 Resp. H335; Aquatic Tox. Chronic 3 H412;
Benzenesulfonic Acid, 4 C10-C13 sec-Alkyl Derivatives	85536-14-7	287-494-3	1.5 to 2	Acute Tox. 4 H302; Skin Corr. 1C H314; Aquatic Tox. Chronic 3 H412

**SECTION 4. First aid measures.**

**4.1 Description of first aid measures**

After inhalation: If inhaled, remove to fresh air. If not breathing, give artificial respiration. If breathing is difficult, give oxygen. Get medical attention immediately.

After skin contact: Flush skin with water. Cover the irritated skin with an emollient. Remove contaminated clothing and shoes. Cold water may be used. Wash clothing before reuse.

After eye contact: Check for and remove any contact lenses. In case of contact, immediately flush eyes with plenty of water for at least 15 minutes. Seek medical attention if irritation develops.

After swallowing: Do NOT induce vomiting unless directed to do so by medical personnel. If large quantities of this material are swallowed, call a physician immediately. Loosen tight clothing such as a collar, tie, belt or waistband. Never give anything by mouth to an unconscious person.

**4.2 Most important symptoms and effects, both acute and delayed**  
Difficultly breathing, Skin irritation, Eye irritation  
Do NOT induce vomiting unless directed to do so by medical personnel. If large quantities of this material are swallowed, call a physician immediately.

**4.3 Indication of any immediate medical attention and special treatment needed.** No special treatment is required

**SECTION 5. Firefighting measures.**

**5.1 Extinguishing media**  
Suitable extinguishing agents: SMALL FIRE: Use DRY chemical powder. LARGE FIRE: Use water spray, fog or foam. Do not use water jet.

**5.2 Special hazards arising from the substance or mixture:** When heated to decomposition it emits toxic fumes of sulfur oxides, and sodium oxide.

**5.3 Advice for firefighters**  
Protective equipment: Wear appropriate PPE for fire conditions including self-contained breathing apparatus for firefighting if necessary. Use water spray to cool unopened containers.

SDS: EB17

Revision nr.2  
Dated 04/21/2019  
Page n. 3 / 6

**ENVIROLOGIX**

**SECTION 6. Accidental release measures.**

**6.1 Personal precautions, protective equipment and emergency procedures:** Use personal protective equipment. Avoid dust formation. Avoid breathing vapors, mist or gas. Assure adequate ventilation. Remove all sources of ignition. Evacuate personnel to a safe area. Avoid breathing dust.

**6.2 Environmental precautions:** Prevent further leakage or spillage if safe to do so. Do not let product enter drains. Discharge into the environment must be avoided.

**6.3 Methods and material for containment and clean up:** Sweep up and shovel. Prevent entry into sewers, dike if needed. Eliminate all ignition sources. Call for assistance on disposal. Finish cleaning by spreading water on the contaminated surface and allow to evacuate through the sanitary system.

**6.4 Reference to other sections:** For safe handling refer to Section 7. For information on PPE refer to Section 8. For disposal, refer to Section 13.

**SECTION 7. Handling and storage.**

**7.1 Precautions for safe handling:** Keep away from heat. Keep away from sources of ignition. Prevent electrostatic buildup. Do not ingest. Do not breathe dust. Wear suitable protective clothing. In case of insufficient ventilation, wear suitable respiratory equipment. If ingested, seek medical advice immediately and show the container or the label. Avoid contact with skin and eyes.

**7.2 Conditions for safe storage, including any incompatibilities:** Keep away from incompatibles such as oxidizing agents. Keep container tightly closed. Keep container in a cool, well-ventilated area.

**7.3 Specific end use(s):** Besides the uses described in Section 1.2 there are no other specific uses

**SECTION 8. Exposure controls/personal protection.**

**8.1 Exposure controls**  
Additional information about design of technical systems: None required

**Exposure limits**  
Components with limit values that require monitoring at the workplace:

Chemical	Exposure Limits
Sodium Lauryl Sulfate	OSHA Observe limits for particulate not otherwise regulated: 15 mg/m <sup>3</sup> total dust, 5 mg/m <sup>3</sup> respirable fraction (OSHA PEL) 10 mg/m <sup>3</sup> inhalable particulate, 3 mg/m <sup>3</sup> respirable particulate. (ACGIH TLV) EH40/2005 Inhalable dust: 10mg/m <sup>3</sup> , Respirable dust: 4mg/m <sup>3</sup>

**Exposure controls - Engineering Controls:** Facilities using or storing this material should be equipped with an eyewash and safety shower. Provide local exhaust or general dilution ventilation.

**Personal protective equipment**  
Breathing equipment: Appropriate respiratory protection should be determined according to local conditions using risk analysis protocols. An approved disposable air-purifying particulate respirator may be used as a backup to engineering controls. Always use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU).

Protection of hands: Handle with gloves. Gloves must be inspected prior to use. Use proper glove removal technique (without touching glove's outer surface) to avoid skin contact with this product. Dispose of contaminated gloves after use in accordance with applicable laws and good laboratory practices. Wash and dry hands.

SDS: EB17

Revision nr.2  
Dated 04/21/2019  
Page n. 5 / 6

**ENVIROLOGIX**

**SECTION 12. Ecological Information.**

**12.1 Toxicity: Sodium Lauryl Sulfate**  
Aquatic toxicity: Note: Aquatic Toxicity of mixture is based on Sodium Lauryl Sulfate.

Aquatic toxicity LC50	Effect dose	Exposure	Species
Acute fish toxicity	10.2-22.8 mg/l	96 hours	Pimephales promelas
Acute daphnia toxicity	1.8 mg/l	48 hours	daphnia magna
Acute algae toxicity	117 mg/l	96 hours	Pseudokirchneriella subcapitata
	53 mg/l	96 hours	Desmodesmus subspicatus
	30-100 mg/l	96 hours	Desmodesmus subspicatus

**12.2 Persistence and degradability:** Biodegradability Result: 90 % - Readily biodegradable. Ratio BOD/ThBOD 95.9 %

**12.3 Bio accumulative potential:** Cyprinus carpio (Carp) - 72 h. Bioconcentration factor (BCF): 3.9 - 5.3

**12.4 Mobility in soil:** Not available

**12.5 Results of PBT and vPvB assessment:** Not available as a chemical safety assessment, not required/not conducted.

**12.6 Other adverse effects** No others listed.

**SECTION 13. Disposal considerations.**

Waste treatment methods/ Uncleaned packaging: Dispose of contents and containers in accordance with local, state and federal regulations.

**SECTION 14. Transport information.**

**14.1 UN-Number DOT, ADR, ADN, IMDG, IATA:** UN2926

**14.2 UN proper shipping name DOT, ADR, ADN, IMDG, IATA:** FLAMMABLE SOLID, TOXIC, ORGANIC, N.O.S. (Sodium dodecyl sulfate)

**14.3 Transport hazard class(es)**  
Class (DOT, ADR, ADN, IMDG, IATA): 4.1 (6.1)

**14.4 Packing group (DOT, ADR, IMDG, IATA):** PG111

**14.5 Environmental hazards**  
Marine pollutant: Not applicable

**14.6 Special precautions for user:** Not applicable.

**14.7 Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code** Not applicable.

SDS: EB17

**ENVIROLOGIX** Revision nr.2 Dated 04/21/2019 Page n. 4 / 6

The selected protective gloves have to satisfy the specifications of EU Directive 89/686/EEC and the standard EN 374 derived from it.

Eye protection: Safety glasses with side shields; goggles. Use equipment for eye protection tested and approved under appropriate government standards such as NIOSH (US) or EN 166(EU). Eye and face protection regulations are described by OSHA (US) in 29 CFR 1910.133. Do not wear contact lenses when working with chemicals.

**SECTION 9. Physical and chemical properties.**

**9.1 Information on basic physical and chemical properties**

Extraction Reagent Powder (EB17) – no CAS number

Appearance: Solid –Powder, White  
 Odor: Odorless  
 Odor threshold: not applicable  
 pH: 9.5 (1% sol/water)  
 Melting point/freezing point: No data available  
 Initial boiling point and boiling range: No data available  
 Flash point: No data available  
 Evaporation rate: No data available  
 Flammability(solid, gas): May be combustible at high temperature  
 Upper/lower flammability or explosive limits: No data available  
 Vapor pressure: No data available  
 Vapor density: No data available  
 Relative density: No data available  
 Solubility(ies): Soluble in water  
 Partition coefficient: n-octanol/water: No data available  
 Auto-Ignition Temperature: No data available  
 Decomposition temperature: No data available  
 Viscosity: No data available  
 Explosive properties: No data available  
 Oxidizing Properties: Not applicable

**9.2 Other information** None

**SECTION 10. Stability and reactivity.**

**10.1 Reactivity:** Not self-reactive.  
**10.2 Chemical stability:** Stable under normal temperatures and pressures  
**10.3 Possibility of hazardous reactions :** Reaction with strong oxidizers may cause fire.  
**10.4 Conditions to avoid :** Heat, flames, and sparks  
**10.5 Incompatible materials:** Oxidizing agents (eg bleach).  
**10.6 Hazardous decomposition products:** Carbon monoxide, carbon dioxide, sulfur oxides, carbon dioxide, nitrogen oxides, silicone Oxides.

**SECTION 11. Toxicological information.**

Acute effects (toxicity tests):

Sodium lauryl sulfate - 151-21-3		
Acute oral toxicity	LD50= 1200 mg/kg	rat
Acute dermal toxicity	LD50= > 2000 mg/kg	rabbit
Acute inhalation toxicity	LC50= 3900 mg/m3, 1hour	rat

Sensitization: No sensitizing effects known

Additional toxicological information: CMR (carcinogenicity, mutagenicity and toxicity for reproduction) – no CMR effects.

SDS: EB17

**ENVIROLOGIX** Revision nr.2 Dated 04/21/2019 Page n. 6 / 6

**SECTION 15. Regulatory information.**

**15.1 Safety, health and environmental regulations**

HMIS Classification (US): Health hazard: 2, Flammability: 1, Physical Hazards: 0  
 NFPA Rating (US): Health hazard: 2, Fire: 1, Reactivity Hazard: 0

US Federal Regulations  
 TSCA: TSCA 8(b) inventory: Sodium lauryl sulfate  
 Health and Safety Reporting List: Listed  
 CERCLA: Listed  
 SARA Section 302 (Extremely Hazardous Substances): Not listed  
 Clean Air Act: Not listed  
 Clean Water Act: Not listed  
 OSHA: Not listed

European International Regulations  
 European labeling in accordance with EC Directives: This product is on the European Inventory of Existing Commercial Chemical Substances (EINECS No. 205-788-1)  
 Listed  
 CLASS D-2B: Material causing other toxic effects (TOXIC).

Canada – DSL/NSL  
 Canada – WHMIS  
 China: Listed on National Inventory. Japan: Listed on National Inventory (ENCs). Korea: Listed on National Inventory (KECI). Philippines: Listed on National Inventory (PICCS). Australia: Listed on AICS.

**15.2 Chemical safety assessment** Not carried out.

**SECTION 16. Other information.**

This information is true based on our present knowledge. However, EnviroLogix makes no representation of its accuracy or completeness. Persons receiving this information must exercise their independent judgment in determining the product's safety and suitability for its intended use. This document shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship EHS Department EnviroLogix Inc.

**Code Definitions:**

H228 Flammable solid.  
 H302 + H322 Harmful if swallowed or inhaled  
 H315 Causes skin irritation.  
 H318 Causes serious eye damage.  
 H335 May cause respiratory irritation.  
 H412 Harmful to aquatic life with long lasting effects.

P264 Wash hands thoroughly after handling.  
 P273 Avoid release to the environment.  
 P300 Wear protective gloves/ eye protection.  
 P301 + P312 IF SWALLOWED: Call a POISON CENTER/doctor/physician if you feel unwell.  
 P304 + P340 IF INHALED: Remove to fresh air and keep comfortable for breathing.  
 P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes; remove contact lenses if present and easy to do. Continue rinsing.  
 P403 + P233 Store in a well ventilated place. Keep container tightly closed.

SDS: EB17

**ENVIROLOGIX** Revision nr.2 Dated 04/29/2019 Page n. 1 / 6

**Safety data sheet**

**SECTION 1. Identification of the substance/mixture and of the company/undertaking**

**1.1 Product identifier**  
 Trade name: DB 5 Dilution Buffer  
 Part number: 11130, 11665, 12495 (KR-266)

**1.2 Relevant identified uses of the substance or mixture and uses advised against application of the substance / the preparation :**  
 Laboratory chemicals; kit component. Not to be used for purposes other than those specified in product literature.

**1.3 Details of the supplier of the safety data sheet**  
 Manufacturer/Supplier: EnviroLogix Inc., 500 Riverside Industrial Pkwy., Portland ME 04103, USA  
 Phone: (207) 797-0300

**1.4 Emergency telephone number:** (207) 797-0300 Technical Service

**SECTION 2. Hazards identification.**

**2.1 Classification of the substance or mixture**  
 Classification according to 29CFR 1910.1200: Eye Damage Category 1  
 Aquatic Toxic, Chronic Category 2

**2.2 Label elements**  
 Labeling according to 29CFR 1910.1200:

Pictogram:

Signal word: Warning

Hazard Statements: H318 Causes serious eye damage  
 H411 Toxic to aquatic life with long lasting effects

Precautionary Statements: P264 Wash hands thoroughly after handling  
 P280 Wear protective gloves/protective clothing/eye protection/face protection  
 P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses if present and easy to do. Continue rinsing.  
 P337+P313 IF eye irritation persists: Get medical attention/advice

**2.3 Other Statements** Restricted to professional users

SDS DB5 Dilution Buffer

**ENVIROLOGIX** Revision nr.2 Dated 04/29/2019 Page n. 2 / 6

**SECTION 3. Composition/information on ingredients.**

Chemical name	CAS No	EC No	Classification According to 29CFR 1910.1200	Amount (%)
Sodium Tetraborate Decahydrate	1303-96-4	215-540-4	H360 Rep 1B	< 3 %
p-tertiary Octylphenoxy polyethyl alcohol (Trilon X-100)	9002-93-1		H302 Acute Tox. Oral 4 H315 Skin Irrit. 2 H318 Eye Dam. 1 H411 Aquatic Chronic 2	1 %
Surfynol	9014-85-1		H315 Skin irritation 2 H318 Eye damage 1 H335 STOT SE 3	2 %
1,2-Benzisothiazolin-3-one (Proxel-GXL)	2634-33-5	220-120-9	H302 Acute Tox. 4; H315 Skin Irrit. 2 H317 Skin Sens. 1 (C <sub>2</sub> : 0.05%) H318 Eye Dam. 1; H400 Aquatic Acute 1	0.048 %

**SECTION 4. First aid measures.**

**4.1 Description of first aid measures**

After inhalation: In case of inhalation. Remove to fresh air. If not breathing give artificial respiration. Get medical attention immediately.  
 After skin contact: In case of skin contact. Remove contaminated clothing and shoes immediately. Wash affected area with mild soap or detergent for at least 10 minutes or until no evidence of chemical remains.  
 After eye contact: In case of eye contact, immediately flush eyes with plenty of water for at least 15 minutes. Lifting eyelids occasionally, until no evidence of chemical remains. Get medical attention immediately.  
 After swallowing: In case of ingestion. DO NOT induce vomiting unless directed to do so by medical personnel. Never give anything by mouth to an unconscious person. Call a physician immediately.

**4.2 Most important symptoms and effects, both acute and delayed:** None

**4.3 Indication of any immediate medical attention and special treatment needed:** None

**SECTION 5. Firefighting measures.**

**5.1 Extinguishing media:** CO2, extinguishing powder or water spray. Fight larger fires with water spray or alcohol resistant foam.

**5.2 Special hazards arising from the substance or mixture:** None

**5.3 Advice for firefighters:** Wear protective gear appropriate for fire conditions including respiratory protective gear.

SDS DB5 Dilution Buffer

Revision nr.2  
Dated 04/29/2019  
Page n. 3 / 6

**SECTION 6. Accidental release measures.**

**6.1 Personal precautions, protective equipment and emergency procedures:** In the case of spilled mixture wear gloves to prevent skin contact. In the case of a large spill, additional protection is recommended.

**6.2 Environmental precautions:** Do not discharge mixture to sewer system or waterways.

**6.3 Methods and material for containment and cleanup:** Absorb in paper towel and discard in appropriate waste. Clean with water afterwards. Large spills may be neutralized with dilute solutions of sodium carbonate or calcium oxide.

**6.4 References to other sections:** For safe handling refer to Section 7. For information on PPE refer to Section 8. For disposal refer to Section 13

**SECTION 7. Handling and storage.**

**7.1 Precautions for safe handling:** Practice good chemical hygiene when handling. Avoid contact with eyes, skin, and clothing.

**7.2 Conditions for safe storage, including any incompatibilities:** Store in tightly closed, non-metal container, in a corrosive compatible area. Prevent direct sunlight and heat. Store in well aired storage rooms.

**7.3 Specific end uses(s):** Apart from the uses mentioned in section 1.2., no other specific uses are stipulated

**SECTION Exposure controls/personal protection.**

**8.1 Exposure limits:**  
Components with limit values that require monitoring at the workplace:

	EH40/2005	OSHA
<b>Sodium Tetraborate Decahydrate</b>	8 Hr TWA = 5mg/m <sup>3</sup>	8 Hr TWA = 10 mg/m <sup>3</sup>

**8.2 Exposure Controls:**

8.2.1 Engineering controls: Facilities using this mixture should be equipped with an eyewash and safety shower. Use general or local exhaust ventilation to keep airborne concentrations below permissible exposure limits.

8.2.2 General protective and hygienic measures:  
The usual precautionary measures should be adhered to when handling chemicals.

**Eye Protection:** Safety glasses with side shields, goggles. Use equipment for eye protection tested and approved under appropriate government standards such as NIOSH (US) or EN 166 (EU). Eye and face protection regulations are described by OSHA (US) in 29CFR1910.133. Do not wear contact lenses when working with chemicals

**Hand Protection:** Handle with gloves. Gloves must be inspected prior to use. Use proper glove removal technique (without touching glove's outer surface) to avoid skin contact with this product. Dispose of contaminated gloves after use in accordance with applicable laws and good laboratory practices. Wash and dry hands. The selected protective gloves have to satisfy the specifications of EU Directive 89/686/EEC and the standard EN 374 derived from it.

**Breathing Equipment:** Appropriate respiratory protection should be determined according to local conditions using risk analysis protocols. An approved disposable air purifying particulate respirator may be used as a backup to engineering controls. Always use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU).

8.2.3 Environmental exposure controls: Contain spills, do not allow into environment

SDS DB5 Dilution Buffer

Revision nr.2  
Dated 04/29/2019  
Page n. 4 / 6

**SECTION 9. Physical and chemical properties.**

**9.1 Information on basic physical and chemical properties:**

a) Appearance: Clear liquid, colorless to slight yellow.  
b) Odor: None  
c) Odor Threshold: No Data Available  
d) pH: 8.6  
e) Melting point/freezing point: No Data Available  
f) Boiling point/Boiling range: No Data Available  
g) Flash point: Not applicable.  
h) Evaporation rate: No Data Available  
i) Flammability (solid, gaseous): No Data Available  
j) Upper/lower flammability or explosive limits: No Data Available  
k) Vapor pressure: No Data Available  
l) Vapor density: No Data Available  
m) Relative density: No Data Available  
n) Solubility(ies): Fully miscible, water.  
o) Partition Coefficient: n-Octanol/water: No Data Available  
p) Auto-ignition temperature: No Data Available  
q) Decomposition temperature: No Data Available  
r) Viscosity: No Data Available  
s) Explosive properties: No Data Available  
t) Oxidizing properties: No Data Available

**9.2 Other information:** No further relevant information available.

**SECTION 10. Stability and reactivity.**

**10.1 Reactivity:** No data available

**10.2 Chemical Stability:** Stable under normal temperatures and pressures.

**10.3 Possibility of hazardous reactions:** Under normal conditions of storage and use, hazardous reactions will not occur.

**10.4 Conditions to avoid:** No specific data

**10.5 Incompatible materials:** No Data Available.

**10.6 Hazardous decomposition products:** Under normal conditions of storage and use, hazardous decompositions products should not be produced.

**SECTION 11. Toxicological information.**

**Information on Toxicological Effects**

**Triton X-100**  
Acute toxicity: Oral LD50 –Rat- 1800mg/kg  
Dermal LD50- Rabbit- 8000 mg/kg

Sensitization: No sensitizing effects known

CMR (carcinogenicity, mutagenicity and toxicity for reproduction) effects: No CMR effects.

Additional toxicological information: No Additional Information

SDS DB5 Dilution Buffer

Revision nr.2  
Dated 04/29/2019  
Page n. 5 / 6

**SECTION 12. Ecological information.**

**12.1 Toxicity:** Fish: LC50 Pimephales promelas (fathead minnow) – 8.9mg/l – 96.0 hr  
Daphnia: EC50 – Daphnia – 26 mg/l – 48 hr

**12.2 Persistence and degradability :** No Data Available

**12.3 Bio accumulative potential:** No Data Available

**12.4 Mobility in soil :** No Data Available

**12.5 Results of PBT and vPvB assessment:** Not available as a chemical safety assessment, not required/not conducted.

**12.6 Other adverse effects:** No Data Available

**SECTION 13. Disposal considerations.**

Waste treatment methods: Contact a licensed professional waste disposal service to dispose of this material. Disposal of surplus or waste solutions must be in accordance with applicable local, state, and national laws and regulations.

**SECTION 14. Transport information.**

**14.1 UN-Number DOT, ADR, ADN, IMDG, IATA :** Not Hazardous for Transport

**14.2 UN proper shipping name DOT, ADR, ADN, IMDG, IATA :** Not Hazardous for Transport

**14.3 Transport hazard class(es) DOT, ADR, ADN, IMDG, IATA):** Not Hazardous for Transport

**14.4 Packing group (DOT, ADR, IMDG, IATA):** Not Hazardous for Transport

**14.5 Environmental hazards** No environmental hazard.

**14.6 Special precautions for user :** None

**14.7 Transport in bulk according to Annex II of MARPOL73/78 and the IBC code:** No information available.

**SECTION 15. Regulatory information.**

**15.1 Safety, health, and environmental regulations**

**US Federal Regulations**  
OSHA SARA 313 Not a hazardous material  
Not listed

**US State Regulations**  
**European/International Regulations**  
European labeling in accordance with EC Directives Not hazardous according to European directives

**15.2 Chemical Safety Assessment** Not carried out

SDS DB5 Dilution Buffer

Revision nr.2  
Dated 04/29/2019  
Page n. 6 / 6

**SECTION 16. Other information.**

*This information is true based on our present knowledge. However, EnviroLogix makes no representation of its accuracy or completeness. Persons receiving this information must exercise their independent judgment in determining the product's safety and suitability for its intended use. This document shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship*

EHS Department  
EnviroLogix Inc.

**Codes:**

**H302** Harmful if swallowed    **H315** Causes skin irritation    **H317** May cause an allergic skin reaction  
**H318** Causes Serious Eye Damage    **H335** May cause respiratory irritation    **H411** Toxic to Aquatic Life with Long Lasting Effects

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