



# DNABLE® Sample Extraction Set 3

## Set Contains:

- MB8 Extraction Buffer
- MB1 Dilution Buffer
- 1.5 mL clear micro-centrifuge tubes (50) for extraction
- 1.5 mL blue micro-centrifuge tubes (50) for dilution

## Materials Not Provided:

- Sterile Modified Buffered Peptone Water (mBPW), 2% (Salmonella assays)
- Oxoid Novel Enrichment (ONE) broth (Listeria assay)
- Precision pipette(s) capable of delivering 20-1000  $\mu$ L
- Pipette tips
- Incubator capable of  $37 \pm 1^\circ\text{C}$
- Dry heat block capable of  $95$  and  $100 \pm 1^\circ\text{C}$ , with insert suitable for 1.5 mL tubes
- Vortexer
- Micro-centrifuge capable of 12,000 x g
- Timer
- Filter bags for culture

Catalog No. ACC-089

Part #11991

## Intended Use

- When used with DNABLE Kits for *Salmonella* (DF-026 and DF-126): This Set provides for extraction and detection of *Salmonella* DNA from a variety of matrices including the whole shrimp matrix. DNABLE *Salmonella* Supplement (Cat. No. XSALMD550) may be required for some matrices. Contact Technical Support for specific recommendations on other matrix types.
- When used with DNABLE Kit for *Listeria monocytogenes* (DF-019): This Set provides for extraction and detection of *Listeria* DNA from a variety of matrices including cheese and stainless steel surface swabs. Oxoid Novel Enrichment (ONE) broth may be required for some matrices. Contact Technical Support for specific recommendations on other matrix types.

## Intended User

DNABLE assays are designed to be simple and user friendly. It is designed for use by personnel with appropriate training in handling human pathogens, and in Microbiology and Molecular Assay techniques. Training specific to the DNABLE assays will be provided by EnviroLogix; contact Technical Service or visit [envirologix.com](http://envirologix.com) for more information.

## How the Kit Works

An aliquot of MB8 buffer is added to a micro-centrifuge tube followed by a sample of mBPW-enriched (*Salmonella*) or ONE broth-enriched (*Listeria*) culture. The sample is concentrated by centrifugation and heated to enable lysis of cells. A centrifugation and dilution step follows.

## Precautions and Notes

DNABLE is a highly sensitive assay. Therefore the following precautions are recommended to reduce the chance of sample contamination:

- Separated work areas are recommended for each of the following:
  - DNABLE culture
  - DNABLE sample preparation
  - DNABLE amplification and detection
- Clean the work station and pipettes with 10% bleach before and after use
- Do not reuse kit disposables
- Change pipette tips in between samples
- Wear gloves and change between handling of samples
- Avoid delays between sample preparation steps and between sample preparation and DNA amplification
- mBPW and/or ONE broth should be equilibrated at  $37^\circ\text{C}$  before use
- Filter bags should be used during enrichment to minimize particulates
- Enriched cultures should be mixed before sampling
- MB8 is stable for 1 year post manufacture when stored refrigerated at  $4-8^\circ\text{C}$ .
- **Safety:** *Salmonella* and *Listeria* are pathogenic and are classified as a Biosafety Level 2 organisms. Personnel should be appropriately trained and should use personal protective equipment. Laboratories should follow appropriate local safety and environmental regulations and guidelines for containment and disposal as described in the Center for Disease Control and Prevention Manual, "Biosafety in Microbiological and Biomedical Laboratories" <http://www.cdc.gov/biosafety/publications/bmbl5/bmbl.pdf>.

## Enrichment

### Notes

- Collect sample according to your facility's sampling plan.
- Composite samples according to procedures defined for your facility.
- For *Salmonella* assay(s): Prepare 2% mBPW using DNable *Salmonella* Media Supplement.
  - Prepare BPW according to manufacturer's instructions and allow cooling to 30°C before adding supplement.
  - Using aseptic technique, add DNable *Salmonella* supplement to BPW at a ratio of 2 to 100 and mix.
  - Once prepared, mBPW may be stored refrigerated at 4-8°C for up to three months, protected from light.
- For *Listeria* assay: ONE broth as per manufacturer's instructions.

Choose protocol based on sample type. Refer to the key listed with each sample type for applicability:  
(S=*Salmonella* DF-026, S+=*Salmonella* Plus, Lm=*Listeria monocytogenes* DF-019)

Whole Shrimp (S, S+)	Cheese (Lm)	Stainless Steel surface (Lm)
<ol style="list-style-type: none"> <li>1. Suspend 25 g of macerated shrimp per 225 mL of sterile mBPW.</li> <li>2. Hand stomach or otherwise homogenize the sample for 2 minutes (<math>\pm 30</math> seconds).</li> <li>3. Incubate 20-22 hours at <math>37 \pm 1^\circ\text{C}</math>.</li> </ol>	<ol style="list-style-type: none"> <li>4. Suspend 25 g samples in 225 mL sterile ONE broth</li> <li>5. Use of filter culture bags is recommended.</li> <li>6. Stomach or otherwise homogenize the sample for 2 minutes.</li> <li>7. Incubate overnight 16-22 hours at <math>37^\circ\text{C}</math>.</li> </ol>	<ol style="list-style-type: none"> <li>1. Pre-moisten collection sponge with ONE broth.</li> <li>2. Swab the surface using horizontal and vertical motions.</li> <li>3. Place the sponge or swab in 225 mL of ONE broth.</li> <li>4. Stomach by hand for 2 minutes.</li> <li>5. Incubate suspended samples overnight 16-22 hours at <math>37^\circ\text{C}</math>.</li> </ol>

Contact Technical Service for additional matrix information and application guides for other matrices not listed herein

## Sample Preparation

<i>Salmonella</i> assays	<i>Listeria</i> assay
<ol style="list-style-type: none"> <li>1. Pre-heat a dry heat block to <math>95^\circ\text{C}</math>. Verify heat block is holding temperature with <math>\pm 1.5^\circ</math> using a simple thermometer.</li> <li>2. Mix the culture before sampling.</li> <li>3. Transfer 1 mL of culture to a clear 1.5 mL micro-centrifuge tube supplied with the set.</li> <li>4. Centrifuge the tube at <math>10,000 \times g</math> for 5 minutes (<math>\pm 30</math> seconds).</li> <li>5. Remove the supernatant using caution to avoid disturbing the pellet. Leave a small volume remaining (<math>\approx 100 \mu\text{L}</math>) if an obvious pellet is not observed.</li> <li>6. Add <math>100 \mu\text{L}</math> of MB8 buffer to the pellet and briefly vortex to suspend the pellet. <i>IMPORTANT: Shake the MB8 vial before pipetting. This is to ensure that the undissolved particulate matter is homogenously distributed in the buffer.</i></li> </ol>	<ol style="list-style-type: none"> <li>1. Pre-heat a dry heat block to <math>100^\circ\text{C}</math>. Verify heat block is holding temperature with <math>\pm 1.5^\circ</math> using a simple thermometer.</li> <li>2. Mix the culture before sampling.</li> <li>3. Transfer 1 mL of enriched sample to a clear 1.5 mL micro-centrifuge tube.</li> <li>4. Centrifuge the tube at <math>12,000 \times g</math> for 3 minutes (<math>\pm 30</math> seconds).</li> <li>5. Remove the supernatant using caution to avoid disturbing the pellet. Leave a small volume remaining (<math>\approx 10 \mu\text{L}</math>) if an obvious pellet is not observed.</li> <li>6. Re-suspend pellet in 1 mL sterile H<sub>2</sub>O</li> <li>7. Centrifuge the tube at <math>12,000 \times g</math> for 1 minute.</li> <li>8. Discard supernatant.</li> <li>9. Add <math>100 \mu\text{L}</math> of MB8 buffer to the pellet and briefly vortex to suspend the pellet.</li> </ol>

- |  |  |
|--|--|
| <ol style="list-style-type: none"> <li>7. Heat sample in heat block at <math>95 \pm 1.5^\circ\text{C}</math> for 5 minutes.</li> <li>8. Centrifuge the tube at <math>10,000 \times g</math> for 5 minutes (<math>\pm 30</math> seconds).</li> <li>9. Place 90 <math>\mu\text{L}</math> of MB1 into a blue 1.5 mL micro-centrifuge tube using a fresh pipette tip.</li> <li>10. Transfer 10 <math>\mu\text{L}</math> of the supernatant from the second centrifugation into the MB1 and mix gently.</li> <li>11. 5 <math>\mu\text{L}</math> of the sample from step 10 will be used in the subsequent DNable reaction.</li> </ol> | <p><i><b>IMPORTANT:</b> Shake MB8 bottle for 5 seconds before use. This is to ensure that the undissolved particulate matter is homogeneously distributed in the buffer.</i></p> <ol style="list-style-type: none"> <li>10. Heat the sample in the heat block at <math>100 \pm 1.5^\circ\text{C}</math> for 10 minutes.</li> <li>11. Vortex briefly (<math>\sim 5</math> seconds) and centrifuge the tube at <math>12,000 \times g</math> for 1 minute.</li> <li>12. Add 20 <math>\mu\text{L}</math> of the supernatant to 80 <math>\mu\text{L}</math> of MB1 using the blue dilution tubes.</li> <li>13. 5 <math>\mu\text{L}</math> of the diluted sample will be used in the DNable reaction.</li> </ol> |
|--|--|

## LIMITED WARRANTY

EnviroLogix Inc. ("EnviroLogix") warrants the products sold hereunder ("the Products") against defects in materials and workmanship when used in accordance with the applicable instructions for a period not to extend beyond a product's printed expiration date. If the Products do not conform to this Limited Warranty and the customer notifies EnviroLogix in writing of such defects during the warranty period, including an offer by the customer to return the Products to EnviroLogix for evaluation, EnviroLogix will repair or replace, at its option, any product or part thereof that proves defective in materials or workmanship within the warranty period.

ENVIROLOGIX MAKES NO OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING BUT NOT LIMITED TO ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. The warranty provided herein and the data, specifications and descriptions of EnviroLogix products appearing in EnviroLogix published catalogues and product literatures are EnviroLogix' sole representations concerning the Products and warranty. No other statements or representations, written or oral, by EnviroLogix' employees, agents or representatives, except written statements signed by a duly authorized officer of EnviroLogix Inc., are authorized; they should not be relied upon by the customer and are not a part of the contract of sale or of this warranty.

EnviroLogix does not warrant against damages or defects arising in shipping or handling, or out of accident or improper or abnormal use of the Products; against defects in products or components not manufactured by EnviroLogix, or against damages resulting from such non-EnviroLogix made products or components. EnviroLogix passes on to customer the warranty it received (if any) from the maker thereof of such non-EnviroLogix made products or components. This warranty also does not apply to Products to which changes or modifications have been made or attempted by persons other than pursuant to written authorization by EnviroLogix.



For Technical Support  
Contact Us At:

**EnviroLogix**  
500 Riverside Industrial Parkway  
Portland, ME 04103-1486 USA  
**Tel: (207) 797-0300**  
**Toll Free: 866-408-4597**  
**Fax: (207) 797-7533**

e-mail: [dnable@envirologix.com](mailto:dnable@envirologix.com)

Website: [www.envirologix.com](http://www.envirologix.com)

**THIS WARRANTY IS EXCLUSIVE.** The sole and exclusive obligation of EnviroLogix shall be to repair or replace the defective Products in the manner and for the period provided above. EnviroLogix shall not have any other obligation with respect to the Products or any part thereof, whether based on contract, tort, strict liability or otherwise. Under no circumstances, whether based on this Limited Warranty or otherwise, shall EnviroLogix be liable for incidental, special, or consequential damages.

This Limited Warranty states the entire obligation of EnviroLogix with respect to the Products. If any part of this Limited Warranty is determined to be void or illegal, the remainder shall remain in full force and effect.

*EnviroLogix, the EnviroLogix logo, DNable, and the DNable logo are trademarks of EnviroLogix Inc.*

© EnviroLogix 2015