

Sample Preparation

Prepare Solutions &
Calibrators



Extract Sample



Microcentrifuge
Sample




EnviroLogix

QuantiTox™ Kit

QuickGuide

Testing Procedure


1 Add Calibrators &
Sample Extract



2 Add Conjugate




3 Mix & Incubate



4 Empty & Wash




Alternate Wash
4 Method




5 Add Substrate



6 Mix & Incubate



7 Add Stop Solution
Read Results



EnviroLogix Technical Support: Toll Free 1-866-408-4597

© EnviroLogix Inc. 2003