


Interpretation Guide

QuickComb for
QuickScan



Lot: 123 - 09

QuickScan
code

C1 C9 RR C3 TF LP 34 3A

Strip #
Strip Color

1 GREEN 2 BLUE 3 RED 4 ORANGE 5 MAROON 6 GOLD 7 GRAY 8 LT. BLUE

Trade names
and proteins

YieldGard Bt11 MON810 Starlink Cry9C Roundup Ready NK603 Event 603 YieldGard Rootworm Event 863 Cry3B Herculex I Cry1F Liberty Link T25 Herculex RW Cry34 Agrisure RW mCry3A MIR604

Detection
levels

0.8% 1/125 0.125% 1/800 0.5% 1/200 0.5% 1/200 0.5% 1/200 0.5% 1/200 0.5% 1/200 1.0% 1/100

