



Sorbent Systems-online presence of IMPAK Corporation
 2460 East 57th St
 Los Angeles, CA 90058
 Voice: 323-277-4700
 Fax: 323-277-9400

Stand-Up Sample Pouches with ZipSeal™



If using bags for mixing ground samples with water, we recommend using sample bags from SorbentSystems.com. You can call the number above or order online:

http://www.sorbentsystems.com/standupbags_order.html

The webpage starts with a size chart, so you can calculate the size bag that is needed:

Stand Up Pouch - Volume/Capacity Chart (all measurements are approximate)

Style	Dimensions	Dry Volume (oz)	Dry Volume (g)	Water Volume (oz)
AA	3.125" x 5" x 1.0"	1 oz	30 g	3 oz
A	4" x 6.41" x 2.25"	2 oz	60 g	5 oz
B	5" x 8" x 3"	4 oz	140 g	16 oz
C	6" x 9.37" x 3.25"	8 oz	250 g	24 oz
D	6.69" x 11" x 3.5"	12 oz	375 g	32 oz
E	7.625" x 11.75" x 4"	16 oz	500 g	48 oz
G	11" x 16.2" x 5.75"	4 lb	1.8 kg	12 lb

- For CORN testing, the size "C" bag works well for an 800 kernel grind (approx. 240g dry volume plus 360 mL [~12 oz] water)
- For SOYBEAN testing without the QuickStix Reader, the size "B" bag works for most 300 bean grinds (approx. 47 g dry volume plus 235 mL [~8 oz] water)
- For SOYBEAN testing with the QuickStix Reader, the size "D" bag is just big enough for most 300 bean grinds (approx. 47 g dry volume plus 940 mL [~32 oz] water)

Scrolling down a bit gets to the section to order standup bags, with the all-clear bags listed first.

Stand Up Pouches with ZipSeal™

- All Dimensions are O.D. (Outside Dimensions)
- 1000 Per Case

Size	Part Number	Dimensions (inches)	Color	1 Case	2 - 4 Cases	5 - 9 Cases	10 - 14 Cases	15 - 24 Cases	25 + Cases
A	AMP02Z	4 x 6.41 x 2.25	Clear	\$109.00	\$86.00	\$79.00	\$73.00	\$68.00	\$63.00
B	BMP03Z	5 x 8 x 3	Clear	\$143.00	\$109.00	\$100.00	\$92.00	\$86.00	\$81.00
C	CMP325Z	6 x 9.37 x 3.25	Clear	\$152.00	\$118.00	\$108.00	\$100.00	\$94.00	\$86.00
D	DMP35Z	6.69 x 11 x 3.5	Clear	\$238.00	\$202.00	\$186.00	\$178.00	\$165.00	\$148.00



Questions? Call EnviroLogix at 1-866-408-4597