

NEW PRODUCT LINE EXCEEDS 2016 GIPSA SPECIFICATIONS



QuickTox for QuickScan has been improved to exceed the new 2016 performance criteria from USDA/GIPSA. EnviroLogix is committed to bringing you tomorrow's solutions today.

FLEX ADVANTAGES OVER THE ORIGINAL QUICKTOX KITS:

- **Accuracy** – test accuracy is improved by eliminating adverse environmental effects
- **Time savings** – the majority of Flex test kits will feature shorter run times and fewer dilution steps
- **Efficiency** – test a wider variety of matrices utilizing the same QuickTox Strips
- **Wider range** – more accurate results across a broader range of test results.

Paired with the proprietary software included in the QuickScan System, QuickTox strips quantify mycotoxins in the ranges most important to you.

QUICKTOX FOR QUICKSCAN STILL OFFERS:

- **Ease of use** – simple procedure, less “shaking,” fewer pipetting steps, and no drying step – read the strips wet!
- **Quantification** – A unique barcoding system containing lot-specific standard curve data ensures accurate and precise mycotoxin measurement while eliminating the need for calibrators or standards.
- **Traceability** – Instantaneous data and image storage onto a host PC (or network) means test results can be emailed, printed or analyzed for management decisions at any time.